

**tedisel**medical

**N270**

**MAINTENANCE MANUAL**



**tediselmedical.com**

# Content

1. Manufacturer .....	4
2. Security information.....	4
2.1. Injury risk warnings .....	4
2.2. Warnings of risk of damage .....	4
2.3. Additional symbols used in the safety instructions .....	5
2.4. Indication of additional information .....	5
2.5. Proper use of oxygen. ....	5
2.5.1. Oxygen explosion .....	5
2.5.2. Fire hazard.....	5
3. Risks.....	6
3.1. Gas explosion .....	6
3.2. Risk of device malfunction .....	6
3.3. Fire risk.....	6
3.4. Danger of electric shock.....	6
4. Symbols used.....	6
5. Product data.....	9
5.1. Storage conditions.....	9
5.2. Operating conditions.....	9
5.3. Service life .....	9
5.4. Purpose of the product .....	10
6. Maintenance .....	10
6.1. Training.....	10
6.2. Removal and fitting of external covers .....	10
6.2.1. Removal of diffusers.....	10
6.2.2. Removal of top and bottom covers.....	11
6.2.3. Removal of central cover .....	12
6.3. Medical gas supply circuits.....	12
6.4. Electrical and voice and data circuits, lighting .....	12
6.5. Replacement of LED strips and drivers in lighting module s.....	13
6.6. Envelopes and structural elements.....	14
6.7. Inspection plan .....	15
7. Cleaning.....	15
8. Waste management.....	16

# N270

Maintenance Manual

9. Regulations .....16

9.1. Team ranking.....16

9.2. Reference standards .....16

9.3. Electromagnetic compatibility. ....17

# N270

## Maintenance Manual

## 1. Manufacturer

Manufacturer: TEDISEL IBÉRICA S.L.

Address: C/ Sant Lluc, 69-81. 08918 - Badalona (Barcelona) SPAIN

Tel. +34 933 992 058

Fax +34 933 984 547

tedisel@tedisel.com

[www.tediselmedical.com](http://www.tediselmedical.com)



## 2. Security information

Important notes in these operating instructions are marked with graphic symbols and signal words.

### 2.1. Injury risk warnings

Signal words such as DANGER, WARNING or CAUTION describe the degree of risk of injury. The different triangular symbols visually emphasise the degree of danger.



WARNING

Refers to a potentially hazardous situation which, if not avoided, could result in death or serious injury.



CAUTION

Refers to a potential hazard which, if not avoided, may result in minor or slight injury.



DANGER

Refers to an immediate danger which, if not avoided, will result in death or serious injury.



Risk of finger entrapment

### 2.2. Warnings of risk of damage

The signal word WARNING describes the degree of risk of material damage. The triangular symbol visually emphasises the degree of danger.



Damage to surfaces: warns of damage to surfaces due to unsuitable cleaning agents and disinfectants.

**NOTICE**

Refers to a potential hazard which, if not avoided, may cause damage to the equipment.

### 2.3. Additional symbols used in the safety instructions



Fire hazard




Explosion hazard: warns of ignition of explosive gas mixtures.



Dangerous voltage: warns about electric shocks that can cause serious injury or death.

### 2.4. Indication of additional information



A NOTE provides additional information and useful tips for safe and efficient use of the device.

### 2.5. Proper use of oxygen.

#### 2.5.1. Oxygen explosion



Oxygen becomes explosive when it comes into contact with oils, greases and lubricants.

Compressed oxygen presents an explosion hazard:

- Make sure that oxygen and gas outlets are free of oil, greasy materials and lubricants!
- Do not use cleaning agents containing oil, grease or lubricants.

#### 2.5.2. fire hazard



**DANGER:** Escaping oxygen is combustible:

- Open fire, red-hot objects and open light are not allowed when working with oxygen!
- Don't smoke!

### 3. Risks

#### 3.1. Gas explosion



Oxygen becomes explosive when it comes into contact with oils, greases and lubricants.

When in contact with oxygen in the air, medical gases may form an explosive or easily flammable gas mixture. The equipment is not suitable for use in environments containing flammable mixtures of anaesthetics with high concentrations of oxygen or nitrous oxide.

If such high concentrations of flammable mixtures of anaesthetics with oxygen or nitrous oxide occur in the environment of the device, there is a risk of ignition under certain conditions.

#### 3.2. Risk of device malfunction



CAUTION: If a device is connected to the equipment and trips the protection mechanism of the corresponding circuit in the health care facility, other devices connected to the equipment will also be de-energised.

#### 3.3. Fire risk



Plug-in connections for the supply of medical gases must not come into contact with oil, grease or flammable liquids.

#### 3.4. Danger of electric shock



Signal cables (network, audio, video, etc.) must be electrically isolated from equipment and the ends of building connections to prevent contact with currents that can cause serious injury or death.

### 4. Symbols used



Applicable part B



Earth (mass)

	Equipotentiality
	Protective earth (ground)
<b>N</b>	Connection point for neutral conductor
	Nurse call button
	Direct lighting
	Indirect lighting
	Operating instructions
<b>MD</b>	Health Product
	Waste electrical equipment
<b>CE</b>	CE symbol
<b>REF</b>	Product code
<b>UDI</b>	Unique identification code



Serial number



Manufacturer



Date of manufacture



Reference to the instruction manual



Damage to surfaces



Fire hazard



Danger of explosion



Dangerous tension



NOTICE

Notice



Risk of finger entrapment



WARNING

Warning

N270

Maintenance Manual



CAUTION

Caution



DANGER

Danger

## 5. Product data

This manual refers to the N270 model. This model is part of the SICA family.

### 5.1. Storage conditions

The individual packaging of this type of product consists of a bubble wrap on the inside and a cardboard box on the outside. Non-stackable packaging.

Under no circumstances should the product be stored with open or damaged packaging. If the product is inspected on receipt and installation is not carried out within 1 day, the product packaging must be resealed.



NOTICE: Failure to follow these instructions may result in damage to the equipment.

Recommended temperature range: -20 °C to 60 °C

Recommended humidity range: 10 % to 75 %.

Atmospheric pressure: 500 hPa to 1,060 hPa

### 5.2. Operating conditions



NOTICE: Failure to follow these instructions may result in damage to the equipment.

Recommended temperature range: -10 °C to 40 °C

Recommended humidity range: 30 % to 75 %.

Atmospheric pressure: 700 hPa to 1,060 hPa

### 5.3. Service life

The useful life of the SICA family of products is determined by the useful life of the medical gas intakes it incorporates, which is 8 years.

## 5.4. Purpose of the product

These systems have three main distinct functions within the hospital:

- Medical gas services
- Electrical, voice and data services
- Lighting
- Nurse call

They consist of a chassis made of aluminium profiles, which integrates the electrical equipment, call, voice and data systems, and installation and channelling of medical gas outlets.

## 6. Maintenance

### 6.1. Training

Personnel performing maintenance must be trained and qualified by the customer. Persons who:

1. have received the training and are duly registered (at those levels where legal provisions make such registration necessary).
2. have been instructed in the maintenance of this device by means of this instruction manual as a basis.
3. are able to assess the tasks they perform on the basis of their own professional experience and training in relevant safety standards and can recognise the potential hazards involved in the work.

### 6.2. Removal and fitting of external covers

For maintenance of the equipment, the outer covers must be removed.

#### 6.2.1. Removal of diffusers

- Remove the light diffusers ① as shown in Fig. 1. Use a flat-bladed screwdriver with care not to scratch the paint on the side covers.




Fig.1 Removal of upper and lower diffusers

- Put the diffusers back on the equipment and press them into position until the clipping sound is heard.

#### 6.2.2. Removal of top and bottom covers

- Remove the diffusers as described in the previous chapter. You now have access to the side of the upper ① and lower ② covers.
- Remove the top covers ① and bottom covers ② as shown in figure 2 and leave them in a safe place.




Fig.2 Removal of top and bottom coverings

- Place the top and bottom covers back on the equipment by bringing them into position and pressing them together until the clipping sound is heard.

### 6.2.3. Removal of central cover

- Remove the front cover ② using the suction cup ① as shown in figure 3 to reveal the pre-installed gas inlets on the device.




Fig.3 Removing the central cover of the gas rail

- Place the central cover on the device and press it into position until the clipping sound is heard.

### 6.3. Medical gas supply circuits

It is recommended that the equipment be disconnected electrically before servicing.

- Remove the front cover protecting the gas pipes.



See section 6.1 of this manual.

- Check the welding points of the pipeline by wetting these areas with a soapy solution (welding test foam).
- Check the condition of the brackets to ensure that they are securely attached to the gas profile.
- Check the condition of the gas intakes, their actuators, etc.

### 6.4. Electrical and voice and data circuits, lighting



It is recommended that the equipment be disconnected electrically before servicing.

- Remove the top and bottom covers protecting the conduits for electrical, voice and data components and lighting.



See section 6.1 of this manual.

NOTA

For equipment in which an electrical element located on the central rail has been requested, the protective cover of this rail must be removed again.

- Carry out a visual inspection before starting the checks.
- Sockets: Check the voltage at each of the equipment's sockets.
- Lighting: Check the switching on/off from the push buttons on the equipment and/or from the call control. If they do not work properly, see point 6.5 of this manual.
- Voice and data: Check each of the mechanisms of the equipment and call control. To be carried out by the centre's IT and communications staff.

### 6.5. Replacement of LED strips and drivers in lighting modules

If the lighting modules of the N270 system are malfunctioning, both the LED strips ② and the controllers ① must be replaced.



Disconnect the equipment electrically before replacement.

- Remove the diffusers as described in paragraph 6.1.1 of this manual. The lighting module shall be exposed.
- Disconnect the quick connector from the LED strip ②.
- Disconnect the power supply of the controller ① from the terminal strip.
- Unscrew the M4 x16 hex screws ④ DIN 933 releasing the tab ③ holding the controller ① and LED strip ②.




Fig.4 Replacing LED strips and controllers

- Attach the LED strip ② and secure it with an M4 x16 hex screw ④ (the one that is not used to secure the tab ④ that holds the driver).
- Fit the new controller ① and secure it with the tab ③ by screwing in the second hexagonal screw ④.

- Connect the power supply of the controller ① back to the terminal strip.
- Connect the power supply quick connector of the newly installed LED strip ②.
- Check that the lighting module is fixed in position.
- Power up the lighting circuit and perform a test run to check that the lighting module switches on and off.



Contact with live parts can cause an electric shock.

Put the covers back in place.

## 6.6. Envelopes and structural elements

- Perform a visual inspection to detect if any items are not properly secured.



In case of suspicion, carry out a physical check of the elements and refasten them properly.

## 6.7. Inspection plan

Inspection plan		
Important information		
<ul style="list-style-type: none"> <li>• Inspection intervals must be observed</li> <li>• The attached system must be checked by Tedisel or a company authorised by Tedisel.</li> </ul> Following the points and intervals described below.		
Visual Inspection (annual)	NOK	OK
• The parts of the bonded system are not deformed and do not show any damage (scratches, cracks, etc.)*	<input type="checkbox"/>	<input type="checkbox"/>
• The system is free of paint defects*	<input type="checkbox"/>	<input type="checkbox"/>
• Covers are available and correctly fitted*	<input type="checkbox"/>	<input type="checkbox"/>
• All rating plates and labels are available and legible	<input type="checkbox"/>	<input type="checkbox"/>
• All electrical and/or voice and data cables must be inspected for damage and replaced if necessary**	<input type="checkbox"/>	<input type="checkbox"/>
• All copper conduit for medical gases and vacuum must be inspected for damage and replaced if necessary**.	<input type="checkbox"/>	<input type="checkbox"/>
Functional check (annual)	NOK	OK
• The system is securely in place;	<input type="checkbox"/>	<input type="checkbox"/>
• All sockets (if any) are fixed and working properly**.	<input type="checkbox"/>	<input type="checkbox"/>
• Tested protective conductor resistance** (EN 62353)	<input type="checkbox"/>	<input type="checkbox"/>
• Gas intakes (if any) were fixedly mounted and are working properly**.	<input type="checkbox"/>	<input type="checkbox"/>
• Inspection of copper piping for medical gases and vacuum (if any) **	<input type="checkbox"/>	<input type="checkbox"/>
• The lighting (if any) is switched on when the switch is operated.	<input type="checkbox"/>	<input type="checkbox"/>
Observations		
Confirmation of the inspection		
The above-mentioned work, including the necessary adjustments and security testing, has been carried out:		
Date	Name (in capitals)	Signature
* Damaged, deformed or missing components should be replaced as a precautionary measure. In this case, please contact the supplier of the suspension system.		
** If one of the above items is found to be non-compliant during the inspection, the attached system must be shut down immediately as a precautionary measure to ensure that avoid further damage to persons and equipment. Notify the system supplier immediately.		

## 7. Cleaning

Perform this operation with slightly moist cleaning instruments to ensure that no liquid enters the equipment. Since no part or component of the system is invasive, sterilisation is not necessary.



Do not use abrasive or very hard cleaning agents that may cause damage to the exterior coatings, such as disinfectants containing sodium hypochlorite, which is highly corrosive to aluminium.



WARNING: Damage to equipment may occur.

**Formaldehyde-free** disinfectants such as Saint Nebul Aid from Proder Pharma are recommended.

Method of application:

1. Dilute 4 pulses of the valve supplied by the manufacturer per 5 litres of water.
2. Spray the compound on the product and let it react for 15 minutes.
3. Remove with water or soap solution with a wrung out cloth.



Switch off the power supply

Contact with live parts can cause an electric shock.

- Always disconnect the device from the main power supply before cleaning and disinfecting it.
- Do not insert objects into the openings of the device.

## 8. Waste management

Applies WEE2012/19 and RoHS directive 2011/65/EU, amendment 2015/863/EU. The equipment has electrical and electronic components, so it cannot be disposed of as organic waste, but as electrical/electronic waste.

## 9. Regulations

### 9.1. Team ranking

According to the new **MDD** regulation **93/42/EEC** concerning medical devices, this product family is classified as:

- Class IIb, by Annex II, excluding section 4, regulation 11.
- Protection level IP20 according to IEC 60529

Equipment intended for continuous operation.

### 9.2. Reference standards

The device complies with the safety requirements of the following standards and directives:

ISO11197: Medical supply units

IEC 60601-1: Medical electrical equipment. General requirements for basic safety and essential performance.

IEC 60601-1-2: Medical electrical equipment. Part 1-2. General requirements for basic safety and essential performance. Collateral standard. Electromagnetic disturbances.

### 9.3. Electromagnetic compatibility.

According to EN 60601-1-2:2015 this equipment is intended for use in the electromagnetic environment specified below. The user of this equipment must satisfy himself that it is being used in such an environment.

Interference emission measurements	Compliance	Comment
HF emissions according to CISPR 11 standard	Group 1	The supply unit uses HF energy exclusively for its internal OPERATION. Therefore, its HF emissions are minimal and interference with devices in its vicinity is unlikely.
HF emissions according to CISPR 11 standard	Class A	The roof supply unit is suitable for use in non-domestic installations and in installations that are directly connected to the PUBLIC SUPPLY NETWORK, which also supplies residential buildings.
Harmonic emissions according to the standard IEC 61000-3-2	Class A	
Emissions of voltage fluctuations/transients in accordance with the standard IEC 61000-3-3	In accordance with	

Interference resistance	Test level according to IEC 60601	Level of compliance	Environment/Guidelines
Static Electric Discharge (ESD) according to IEC 61000-4-2	±8 kV contact discharge 15 kV aerial discharge	±8 kV contact discharge 15 kV aerial discharge	Floors should be made of wood, concrete or ceramics. If the floor is covered with a synthetic material, the relative air humidity should be at least 30%.
Fast transient electrical interference amplitudes / bursts according to the norm IEC 61000-4-4	±2 kV for power supply cables ±1kV for input and output cables	±2 kV for power supply cables ±1 kV for incoming and outgoing cables	The quality of the supply voltage should be typical for a commercial or hospital environment.

Overvoltage (wave) according to IEC 61000-4- 5	±1 kV phase-to-phase voltage ±2 kV phase to ground voltage	±1 kV phase-to-phase voltage ±2 kV phase to ground voltage	The quality of the supply voltage should be typical for a commercial or hospital environment.
Voltage dips and fluctuations of the supply voltage according to the standard IEC 61000-4- 11	100% of $U_N$ drop for 0.5 period 100% of $U_N$ drop for 1 period 30% of $U_N$ drop for 25 periods  Remark: $U_N$ is the AC mains voltage before applying the test level.	100% $U_N$ drop for 0.5 period 100% of $U_N$ drop for 1 period 30% of $U_N$ drop for 25 periods	The quality of the supply voltage should be typical for a commercial or hospital environment.  If the user of the roof supply unit requires continuous operation even in case of power supply interruptions, it is recommended to supply the roof supply unit from a device with an uninterruptible power supply or a battery.
Short interruptions of the supply voltage according to the standard IEC 61000-4- 11	100% for 5 s  Remark: $U_N$ is the AC mains voltage before applying the test level.		The quality of the supply voltage should be typical for a commercial or hospital environment.  If the user of the roof supply unit requires continuous operation even in case of power supply interruptions, it is recommended to supply the roof supply unit from a device with an uninterruptible power supply or a battery.
Magnetic field for power supply frequencies (50/60 Hz) according to the standard IEC 61000-4-8	30 A/m	30 A/m	The magnetic fields created by the mains frequency should be those of a commercial or hospital environment.

Interference resistance	Level of verification according to IEC 60601	Level of compliance	Environment/Guidelines																																																		
HF interference induced by IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM band	3 Vrms 6 Vrms	AM 1KHz modulation Depth 80% Depth 80% Depth 80% Depth																																																		
HF interference induced by IEC 61000-4-3	<table border="1"> <thead> <tr> <th>RANGE</th> <th>FREQUENCY</th> <th>MODULATION</th> <th>STEP</th> <th>LEVEL</th> </tr> </thead> <tbody> <tr> <td>A</td> <td>80-1000MHz</td> <td>AM 1 kHz Prof: 80%</td> <td>LOG 1%</td> <td>10 V/m</td> </tr> <tr> <td>B</td> <td>1000-2000MHz</td> <td>AM 1 kHz Prof: 80%</td> <td>LOG 1%</td> <td>10 V/m</td> </tr> <tr> <td>C</td> <td>2000-2700MHz</td> <td>AM 1 kHz Prof: 80%</td> <td>LOG 1%</td> <td>10 V/m</td> </tr> <tr> <td>D</td> <td>385MHz</td> <td>PM 18 Hz Cycle: 50%</td> <td>-</td> <td>27 V/m</td> </tr> <tr> <td>E</td> <td>450MHz</td> <td>FM 1 kHz Desv:± 5 kHz</td> <td>-</td> <td>28 V/m</td> </tr> <tr> <td>F</td> <td>810-930MHz</td> <td>PM 18 Hz Cycle: 50%</td> <td>-</td> <td>28 V/m</td> </tr> <tr> <td>G</td> <td>1720-1970MHz</td> <td>PM 217 Hz Cycle: 50%</td> <td>-</td> <td>28 V/m</td> </tr> <tr> <td>H</td> <td>2450MHz</td> <td>PM 217 Hz Cycle: 50%</td> <td>-</td> <td>28 V/m</td> </tr> <tr> <td>I</td> <td>5240-5785MHz</td> <td>PM 217 Hz Cycle: 50%</td> <td>-</td> <td>9 V/m</td> </tr> </tbody> </table>			RANGE	FREQUENCY	MODULATION	STEP	LEVEL	A	80-1000MHz	AM 1 kHz Prof: 80%	LOG 1%	10 V/m	B	1000-2000MHz	AM 1 kHz Prof: 80%	LOG 1%	10 V/m	C	2000-2700MHz	AM 1 kHz Prof: 80%	LOG 1%	10 V/m	D	385MHz	PM 18 Hz Cycle: 50%	-	27 V/m	E	450MHz	FM 1 kHz Desv:± 5 kHz	-	28 V/m	F	810-930MHz	PM 18 Hz Cycle: 50%	-	28 V/m	G	1720-1970MHz	PM 217 Hz Cycle: 50%	-	28 V/m	H	2450MHz	PM 217 Hz Cycle: 50%	-	28 V/m	I	5240-5785MHz	PM 217 Hz Cycle: 50%	-	9 V/m
RANGE	FREQUENCY	MODULATION	STEP	LEVEL																																																	
A	80-1000MHz	AM 1 kHz Prof: 80%	LOG 1%	10 V/m																																																	
B	1000-2000MHz	AM 1 kHz Prof: 80%	LOG 1%	10 V/m																																																	
C	2000-2700MHz	AM 1 kHz Prof: 80%	LOG 1%	10 V/m																																																	
D	385MHz	PM 18 Hz Cycle: 50%	-	27 V/m																																																	
E	450MHz	FM 1 kHz Desv:± 5 kHz	-	28 V/m																																																	
F	810-930MHz	PM 18 Hz Cycle: 50%	-	28 V/m																																																	
G	1720-1970MHz	PM 217 Hz Cycle: 50%	-	28 V/m																																																	
H	2450MHz	PM 217 Hz Cycle: 50%	-	28 V/m																																																	
I	5240-5785MHz	PM 217 Hz Cycle: 50%	-	9 V/m																																																	

Transmitter power rating	Safety distance depending on emission frequency Environment/Guidelines		
	150 kHz to 80 MHz $D = 1,2 P$	80 MHz up to 800 MHz $D = 1,2 P$	800 MHz up to 2.5 GHz $D = 2, 3 P$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

N270

Maintenance Manual