


**tedisel**medical



ANTEA

INTEGRAL DESIGN BED HEAD UNIT

## INTEGRAL DESIGN BED HEAD UNIT



Wall Recessed



Horizontal or Vertical orientation



Colors, HPL or Vinyl upon request



Easy and fast maintenance



Length and specifications upon request



Flush-mounted elements for easy cleaning




Anodized aluminum



Compact model



Durability



### Bringing efficiency to small spaces

**ANTEA** is our fully recesses model designed for hospitalization areas. It is easy to integrate it on HPL panels.

### Some references



Clínica Corachan  
Hospital St. Juan Brugge  
Hospital la Rochelle  
Hospital Costa del Sol  
Gargash Hospital  
Complejo Hospitalario Virgen de la Macarena

Spain  
Belgium  
France  
Spain  
UAE  
Spain



## ANTEA PARTS

**ANTEA** has a front clipped by means of ball retainers and a front panel removal by suction cups.



## QUALITY

Tedisel Management System and all our products are certified with quality seal ISO 9001, ISO 13485 and CE marking granted by notified body TÜV RHEINLAND LGA Products GmbH. Tedisel, as medical device manufacturer, has the medical device manufacturer license number 6205-PS, granted by the Spanish Agency for Medicines and Health Products.

All medical devices developed and manufactured in Tedisel, are regulated by the European Directive 93/42/EEC. All equipment has accreditation and certification of electrical safety and electromagnetic compatibility standards ISO 11197 and EN 60601.

Tedisel Iberica S. L.  
Sant LLuc, 69-81  
08918 Badalona, BARCELONA, Spain  
Tel. +34 93 399 20 58  
info@tedisel.com  
tediselmedical.com



INF-306-EN 09/2023