

# ANTEA

INTEGRAL DESIGN BED HEAD UNIT

## INTEGRAL DESIGN BED HEAD UNIT



Wall Recessed



Horizontal or Vertical  
orientation



Colors, HPL or Vinyl  
upon request



Easy and fast  
maintenance



Length and specifications  
upon request



Flush-mounted elements  
for easy cleaning




Anodized aluminum



Compact model



Durability



### Bringing efficiency to small spaces

**ANTEA** is our fully recesses model designed for hospitalization areas. It is easy to integrate it on HPL panels.


### Some references

Clínica Corachan	Spain
Hospital St. Juan Brugge	Belgium
Hospital la Rochelle	France
Hospital Costa del Sol	Spain
Gargash Hospital	UAE
Complejo Hospitalario Virgen de la Macarena	Spain



## ANTEA PARTS

**ANTEA** has a front clipped by means of ball retainers and a front panel removal by suction cups.



## QUALITY

Tedisel Management System and all our products are certified with quality seal ISO 13485 and CE marking granted by notified body TÜV RHEINLAND LGA Products GmbH. Tedisel, as medical device manufacturer, has the medical device manufacturer license number 6205-PS, granted by the Spanish Agency for Medicines and Health Products.

All medical devices developed and manufactured in Tedisel, are regulated by the European Directive 93/42/EEC. All equipment has accreditation and certification of electrical safety and electromagnetic compatibility standards ISO 11197 and EN 60601.

Tedisel Iberica S. L.  
Sant LLuc, 69-81  
08918 Badalona, BARCELONA, Spain  
Tel. +34 93 399 20 58  
info@tedisel.com  
tediselmedical.com



INF-306-EN 05/2024