

tediselmedical
AN EXPERT INSIGHT INTO HOSPITAL EQUIPMENT



ATLAS

CEILING SUSPENDED BRIDGE

CEILING SUSPENDED BRIDGE



Possibility to benefit from the front and rear sides



Possibility to have sliding vertical trolleys with a large variety of accessories (up to two trolleys per bed)



Rounded sides and geometric lines for optimal cleaning



Colors, HPL, Vinyl, and anodization upon request



Antibacterial paint upon request



Length and specifications upon request



LED lighting
(Direct, Indirect and Night Light)


Effectiveness in Critical Areas

ATLAS is our ceiling suspended bridge designed to optimize space in critical areas and to cover all medical needs.


ATLAS offers different solutions and highlights the capacity to install supply units in its central profile and support for medical material in their trolleys.

- 1 Suspension pipes
- 2 Suspended bridge
- 3 Infusion trolley
- 4 Monitoring trolley
- 5 Shelf for trolley
- 6 Shelf with drawer for trolley


3




4




1




2



5



6



Some references

Hospital Universitario de Jaén	Spain	Vito Fazzi In Leece	Italy
Clínica Fátima, Sevilla	Spain	Burjeel Hospital	EUA
Clínica Iberoamericana, Barranquilla	Colombia	Southern Philippines Medical Center (Davao)	Philippines
Clínica Internacional San Borja, Lima	Peru	Clínica Girona	Spain
Lybian European Hospital	Lybia	Hospital Quirón Marbella	Spain
Centre d'Oncologie Ryad	Morocco		




ACCESSORIES

ATLAS offers different solutions, its main feature is the ability to place various accessories on an optional horizontal rail and vertical trolleys.



Wall Technical rails



IV Poles & Infusion pumps



Shelves & Drawers



Monitor Support Arms



Examination lamps



Medical baskets and more



Mobile and fixed vertical trolleys (Wet and Dry)

QUALITY

Tedisel Management System and all our products are certified with quality seal ISO 13485 and CE marking granted by notified body TÜV RHEINLAND LGA Products GmbH. Tedisel, as a medical device manufacturer, has the medical device manufacturer license number 6205-PS, granted by the Spanish Agency for Medicines and Health Products. All medical devices

developed and manufactured in Tedisel, are regulated by the European Directive 93/42/EEC. All equipment has accreditation and certification of electrical safety and electromagnetic compatibility standards ISO 11197 and EN 60601.

Tedisel Iberica S. L.
Sant Lluç, 69-81
08918 Badalona, BARCELONA, Spain
Tel. +34 93 399 20 58
info@tedisel.com
tediselmedical.com



INF-311-EN

05/2024